

Efficiency in clinics



Up to 50 Hygosuc devices can be connected to the CDS60 central disinfectant station for use in large dental practices and clinics. The CDS60 facilitates feeding of the Orotol® plus concentrate to each Hygosuc for fast and cost-efficient disinfection of the suction system.

Monitoring with VistaSoft Monitor

VistaSoft Monitor software

- Monitoring of the device status via network
- Logging of disinfection cycles
- Increased operational safety and reliability
- Fast help for servicing



Technical data at a glance

Hygosuc	
Voltage (V AC/DC)	24
Current consumption (mA)	Max. 200
Power consumption (VA)	Max. 5
Dimensions of add-on device (H x W x D cm)	32 x 19 x 16
Weight (kg)	5
Water pressure (bar)	4 ± 1
Air pressure (bar)	4 ± 1
Outflow diameter (mm)	Min. 15 mm

DÜRR DENTAL SE
 Höpfigheimer Str. 17
 74321 Bietigheim-Bissingen
 Germany
www.duerrdental.com
info@duerrdental.com



Precise dispensing and disinfection – the Hygosuc from Dürr Dental

Cleaning and disinfection of the suction system made easy.



COMPRESSED AIR
 SUCTION
 IMAGING
 DENTAL CARE
 HYGIENE



Functions of the Hygosuc

The **Hygosuc from Dürr Dental** enables precise, automated cleaning and disinfection of the suction system combined with easy handling, generating substantial time savings for your practice staff. It eliminates the need for routine manual cleaning and disinfection of the suction system. The compact design of the add-on device enables quick and easy retrofitting of the treatment unit. The cleaning and disinfection cycles can be logged via the VistaSoft Monitor software.



Automated cleaning and disinfection

One suction point each for the small and large suction hoses for time-efficient and easy cleaning and disinfection of the suction system.



Integrated flush-cleaning

The recommended suction of min. 200 ml water after each treatment is carried out via the integrated flush-cleaning.

- Key features:**
- Simple, time-efficient handling
 - Proper cleaning and disinfection with precise dispensing
 - Logging of disinfection cycles via VistaSoft Monitor software
 - Low space requirements due to compact design
 - 3 functions in 1 device
 - Flush-cleaning
 - Automated disinfection with Orotol® plus
 - Integrated rinsing unit
 - Made in Germany



The Hygosuc add-on device (also available as a built-in version).

Hygosuc at a glance



- Suction points**
There are several suction points with diameters of 6, 11 and 16 mm.
- Large suction point**
Suction of spray mist
- Small suction point**
Saliva ejector
- Blue suction point**
Flush-cleaning

Button
The two-colour button displays the operating status of the system and indicates when the Orotol® plus disinfection bottle needs changing.

Basic unit
The basic unit is responsible for precise dispensing: Orotol® plus is mixed with water. As an option, the Hygosuc basic unit can also serve as a rinsing unit. Water is continuously supplied to the system during suction. This prevents blood from coagulating and stops deposits that may be difficult to remove from building up in the suction system. The system remains clean and functions reliably.



1 l bottle of Orotol® plus for the Hygosuc add-on device.

Orotol® plus disinfection of suction systems

Foam-free liquid concentrate for simultaneous disinfection, deodorisation, cleaning and care of all suction systems and amalgam separators.

- Key features:**
- For suction systems, recommended by leading manufacturers
 - No development of resistance to germs
 - Certified material compatibility
 - Preserves the value of your suction system
 - Dissolves and disinfects bio-film and prevents the formation of sludge (blood, proteins)
 - Long-term effect
 - Broad range of efficacy